# Package 'ddPCRclust'

November 15, 2024

```
Title Clustering algorithm for ddPCR data
```

**Version** 1.26.0

**Description** The ddPCRclust algorithm can automatically quantify the CPDs of non-orthogonal ddPCR reactions

with up to four targets. In order to determine the correct droplet count for each target, it is crucial to both identify all clusters and label them correctly based on their position. For more information on what data can be analyzed and how a template needs to be formatted, please check the vignette.

```
LazyData true
```

**Depends** R (>= 3.5)

**Imports** plotrix, clue, parallel, ggplot2, openxlsx, R.utils, flowCore, flowDensity (>= 1.13.3), SamSPECTRAL, flowPeaks

Suggests BiocStyle

Collate 'cluster functions.R' 'functions.R' 'ddPCRclust.R'

License Artistic-2.0

URL https://github.com/bgbrink/ddPCRclust

BugReports https://github.com/bgbrink/ddPCRclust/issues

biocViews ddPCR, Clustering

RoxygenNote 6.0.1

 $\textbf{git\_url} \ \ \textbf{https://git.bioconductor.org/packages/ddPCRclust}$ 

git\_branch RELEASE\_3\_20

git\_last\_commit 5098ee0

git\_last\_commit\_date 2024-10-29

**Repository** Bioconductor 3.20

**Date/Publication** 2024-11-14

**Author** Benedikt G. Brink [aut, cre], Justin Meskas [ctb],

Ryan R. Brinkman [ctb]

Maintainer Benedikt G. Brink <br/> <br/>bbrink@cebitec.uni-bielefeld.de>

2 calculateCPDs

# **Contents**

calc	ulateCPDs	Calculates the copies per droplet	
Index			15
	shearCorrection		. 13
	runDensity		. 10
	readTemplate		. 9
	readFiles		. 9
	exportToExcel		. 8
	exportToCSV		. 7
	exportPlots		. 6
	createEnsemble		. 3
	calculateCPDs		. 2

**Description** 

This function takes the results of the clustering and calculates the actual counts per target, as well as the counts per droplet (CPD) for each marker.

#### Usage

```
calculateCPDs(results, template = NULL, constantControl = NULL)
```

# **Arguments**

results The result of the ddPCRclust algorithm.

template The parsed dataframe containing the template.

constantControl

The constant refrence control, which should be present in each reaction. It is used to normalize the data.

# Value

A list of lists, containing the counts for empty droplets, each marker with both total droplet count and CPD, and total number of droplets, for each element of the input list respectively.

```
# Read files
exampleFiles <- list.files(paste0(find.package('ddPCRclust'), '/extdata'), full.names = TRUE)
files <- readFiles(exampleFiles[3])
# To read all example files uncomment the following line
# files <- readFiles(exampleFiles[1:8])
# Read template
template <- readTemplate(exampleFiles[9])</pre>
```

createEnsemble 3

```
# Run ddPCRclust
result <- ddPCRclust(files, template)

# Calculate the CPDs
markerCPDs <- calculateCPDs(result, template$template)</pre>
```

createEnsemble

Create a cluster ensemble

# Description

This function takes the three (or less) clustering approaches of the ddPCRclust package and combines them to one cluster ensemble. See cl\_medoid for more information.

# Usage

```
createEnsemble(dens = NULL, sam = NULL, peaks = NULL, file)
```

# **Arguments**

dens	The result of the flowDensity algorithm as a CLUE partition.
sam	The result of the samSPECTRAL algorithm as a CLUE partition.
peaks	The result of the flowPeaks algorithm as a CLUE partition.
file	The input data. More specifically, a data frame with two dimensions, each dimension representing the intensity for one color.

#### Value

data	The original input data minus the removed events (for plotting)
confidence	The agreement between the different clustering results in percent. If all algorithms calculated the same result, the clustering is likely to be correct, thus the confidence is high.
counts	The droplet count for each cluster.

```
exampleFiles <- list.files(paste0(find.package('ddPCRclust'), '/extdata'), full.names = TRUE)
file <- read.csv(exampleFiles[3])
densResult <- runDensity(file = file, numOfMarkers = 4)
samResult <- runSam(file = file, numOfMarkers = 4)
peaksResult <- runPeaks(file = file, numOfMarkers = 4)
superResult <- createEnsemble(densResult, samResult, peaksResult, file)</pre>
```

4 ddPCRclust

ddPCRclust	ddPCRclust A package for automated quantification of multiplexed ddPCR data

#### **Description**

The ddPCR clust algorithm can automatically quantify the events of ddPCR reaction with up to four markers. In order to determine the correct droplet count for each marker, it is crucial to both identify all clusters and label them correctly based on their position. For more information on what data can be analyzed and how a template needs to be formatted, please check the project repository on github.

This is the main function of this package. It automatically runs the ddPCRclust algorithm on one or multiple csv files containing the raw data from a ddPCR run with up to 4 markers.

# Usage

```
ddPCRclust(files, template, numOfMarkers = 4, sensitivity = 1,
    similarityParam = 0.95, distanceParam = 0.2, fast = FALSE,
    multithread = FALSE)
```

# **Arguments**

files	The input data obtained from the csv files. For more information, please see	
	and the state of t	

readFiles.

template A data frame containing information about the individual ddPCR runs. An ex-

ample template is provided with this package. For more information, please see

readTemplate.

numOfMarkers The number of primary clusters that are expected according the experiment set

up. Can be ignored if a template is provided. Else, a vector with length equal to length(files) should be provided, containing the number of markers used for

the respective reaction.

sensitivity A number between 0.1 and 2 determining sensitivity of the initial clustering,

e.g. the number of clusters. A higher value means the data is divided into more clusters, a lower value means more clusters are merged. This allows fine tuning

of the algorithm for exceptionally low or high CPDs.

similarityParam

If the distance of a droplet between two or more clusters is very similar, it will not be counted for either. The standard it 0.95, i.e. at least 95% similarity. A sensible value lies between 0 and 1, where 0 means none of the 'rain' droplets

will be counted and 1 means all droplets will be counted.

distanceParam When assigning rain between two clusters, typically the bottom 20% are as-

signed to the lower cluster and the remaining 80% to the higher cluster. This parameter changes the ratio, i.e. a value of 0.1 would assign only 10% to the

lower cluster.

fast Run a simpler version of the algorithm that is about 10x faster. For clean data,

this might already deliver very good results. However, is is mostly intended to

get a quick overview over the data.

multithread Distribute the algorithm amongst all CPU cores to speed up the computation.

ddPCRclust 5

#### Value

results

The results of the ddPCRclust algorithm. It contains three fields: data The original input data minus the removed events (for plotting) confidence The agreement between the different clustering results in percent If all parts of the algorithm calculated the same result, the clustering is likely to be correct, thus the confidence is high counts The droplet count for each cluster

#### Usage

The main function of the package is ddPCRclust. This function runs the algorithm with one or multiple files, automatically distributing them amongst all cpu cores using the parallel package (parallelization does not work on windows). Afterwards, the results can be exported in different ways, using exportPlots, exportToExcel and exportToCSV. Once the clustering is finished, copies per droplet (CPD) for each marker can be calculated using calculateCPDs.

These functions provide access to all functionalities of the ddPCRclust package. However, expert users can directly call some internal functions of the algorithm, if they find it necessary. Here is a list of all available supplemental functions:

```
runDensity
runSam
runPeaks
createEnsemble
```

#### Author(s)

 $\textbf{Maintainer}: Benedikt \ G. \ Brink \ \verb|\color| birink@cebitec.uni-bielefeld.de>$ 

Other contributors:

- Justin Meskas < jmeskas@bccrc.ca> [contributor]
- Ryan R. Brinkman <rbrinkman@bccrc.ca> [contributor]

#### See Also

Useful links:

- https://github.com/bgbrink/ddPCRclust
- Report bugs at https://github.com/bgbrink/ddPCRclust/issues

```
# Read files
exampleFiles <- list.files(paste0(find.package('ddPCRclust'), '/extdata'), full.names = TRUE)
files <- readFiles(exampleFiles[3])
# To read all example files uncomment the following line
# files <- readFiles(exampleFiles[1:8])
# Read template
template <- readTemplate(exampleFiles[9])
# Run ddPCRclust
result <- ddPCRclust(files, template)
# Plot the results</pre>
```

6 exportPlots

```
library(ggplot2)
p <- ggplot(data = result$B01$data, mapping = aes(x = Ch2.Amplitude, y = Ch1.Amplitude))
p <- p + geom_point(aes(color = factor(Cluster)), size = .5, na.rm = TRUE) +
    ggtitle('B01 example')+theme_bw() + theme(legend.position='none')
p</pre>
```

exportPlots

Plot the algorithms results with ggplot2

#### **Description**

A convinience function that takes the results of the ddPCRclust algorithm, plots them using the ggplot2 library and a custom colour palette and saves the plots to a folder.

#### Usage

```
exportPlots(data, directory, annotations, format = "png", invert = FALSE)
```

#### **Arguments**

data

The result of the ddPCRclust algorithm

directory

The parent directory where the files should saved. A new folder with the experiment name will be created (see below).

Some basic metadata about the ddPCR reaction. If you provided ddPCRclust a template, this paramater can be filled with the corresponding field in the result. Otherwise, you have to provide a character vector containing a name and the the color channels, e.g. c(Name='ddPCR\_01-04-2017', Ch1='HEX', Ch2='FAM')

format

Which file format to use. Can be either be a device function (e.g. png), or one of 'eps', 'ps', 'tex' (pictex), 'pdf', 'jpeg', 'tiff', 'png', 'bmp', 'svg' or 'wmf' (windows only). See also ggsave

Invert the axis, e.g. x = Ch2. Amplitude, y = Ch1. Amplitude

#### Value

None

invert

```
# Read files
exampleFiles <- list.files(paste0(find.package('ddPCRclust'), '/extdata'), full.names = TRUE)
files <- readFiles(exampleFiles[3])
# To read all example files uncomment the following line
# files <- readFiles(exampleFiles[1:8])
# Read template
template <- readTemplate(exampleFiles[9])
# Run ddPCRclust
result <- ddPCRclust(files, template)
# Export the plots</pre>
```

exportToCSV 7

```
dir.create('./Results')
exportPlots(data = result, directory = './Results/', annotations = result$annotations)
```

exportToCSV

Export the algorithms results to a csv file

# **Description**

A convinience function that takes the results of the droplClust algorithm and exports them to a csv file.

#### Usage

```
exportToCSV(data, directory, annotations, raw = FALSE)
```

#### **Arguments**

data The result of the ddPCRclust algorithm

directory The parent directory where the files should saved. A new folder with the exper-

iment name will be created (see below).

annotations Some basic metadata about the ddPCR reaction. If you provided ddPCRclust a

template, this paramater can be filled with the corresponding field in the result. Otherwise, you have to provide a character vector containing a name and the the color channels, e.g. c(Name='ddPCR\_01-04-2017', Ch1='HEX', Ch2='FAM')

raw Boolean which determines if the annotated raw data should be exported along

with the final counts. Basically, a third column will be added to the original data, which contains the cluster number to which this point was assigned to. Useful for example to visualize the clustering later on. (Warning: this can take a while!)

#### Value

None

```
# Read files
exampleFiles <- list.files(paste0(find.package('ddPCRclust'), '/extdata'), full.names = TRUE)
files <- readFiles(exampleFiles[3])
# To read all example files uncomment the following line
# files <- readFiles(exampleFiles[1:8])

# Read template
template <- readTemplate(exampleFiles[9])

# Run ddPCRclust
result <- ddPCRclust(files, template)

# Export the results
dir.create('./Results')
exportToCSV(data = result, directory = './Results/', annotations = result$annotations)</pre>
```

8 exportToExcel

			-
exp	ort	ToExc	:e1

Export the algorithms results to an Excel file

#### **Description**

A convinience function that takes the results of the droplClust algorithm and exports them to an Excel file.

#### Usage

```
exportToExcel(data, directory, annotations, raw = FALSE)
```

# **Arguments**

data The result of the ddPCRclust algorithm

directory The parent directory where the files should saved. A new folder with the exper-

iment name will be created (see below).

annotations Some basic metadata about the ddPCR reaction. If you provided ddPCRclust a

template, this paramater can be filled with the corresponding field in the result. Otherwise, you have to provide a character vector containing a name and the the color channels, e.g. c(Name='ddPCR\_01-04-2017', Ch1='HEX', Ch2='FAM')

raw Boolean which determines if the annotated raw data should be exported along

with the final counts. Basically, a third column will be added to the original data, which contains the cluster number to which this point was assigned to. Useful for example to visualize the clustering later on. (Warning: this can take a while!)

# Value

None

```
# Read files
exampleFiles <- list.files(paste0(find.package('ddPCRclust'), '/extdata'), full.names = TRUE)
files <- readFiles(exampleFiles[3])
# To read all example files uncomment the following line
# files <- readFiles(exampleFiles[1:8])

# Read template
template <- readTemplate(exampleFiles[9])

# Run ddPCRclust
result <- ddPCRclust(files, template)

# Export the results
dir.create('./Results')
exportToExcel(data = result, directory = './Results/', annotations = result$annotations)</pre>
```

readFiles 9

readFiles	Read the csv files from your disk
-----------	-----------------------------------

# **Description**

This function reads the raw csv files for ddPCRclust from disk and returns the experiment data. Please refer to the vignette for more information on how these files need to be formatted.

#### Usage

```
readFiles(files)
```

#### **Arguments**

files The input file(s), specifically csv files. Each file represents a two-dimensional

data frame. Each row within the data frame represents a single droplet, each

column the respective intensities per colour channel.

#### Value

files A data frame composed of the experiment data

ids The file ids, e.g. A01, A02, etc.

# **Examples**

```
# Read files
exampleFiles <- list.files(paste0(find.package('ddPCRclust'), '/extdata'), full.names = TRUE)
files <- readFiles(exampleFiles[1:8])</pre>
```

readTemplate

Read a template file from disk

# Description

This function reads a template file for ddPCRclust from disk and returns a run template and annotations. Please refer to the vignette for information on how this file need to be formatted.

# Usage

```
readTemplate(template)
```

# **Arguments**

template

A csv file containing information about the individual ddPCR runs. An example template is provided with this package. For more information, please check the vignette or the repository on github.

10 runDensity

#### Value

annotations The metatdata provided in the header of the template. It contains four fields:

Name The name given to this ddPCR experiment

Ch1 Color channel 1 (usually HEX) Ch2 Color channel 2 (usually FAM)

descriptions Additional descriptions about this ddPCR experiment (e.g. date,

exprimentor, etc.)

template A parsed dataframe containing the template.

# **Examples**

```
# Read template
exampleFiles <- list.files(paste0(find.package('ddPCRclust'), '/extdata'), full.names = TRUE)
template <- readTemplate(exampleFiles[9])</pre>
```

runDensity

Find the clusters using flowDensity

#### **Description**

Use the local density function of the flowDensity package to find the cluster centres of the ddPCR reaction. Clusters are then labelled based on their rotated position and lastly the rain is assigned.

#### Usage

```
runDensity(file, sensitivity = 1, numOfMarkers, missingClusters = NULL,
  similarityParam = 0.95, distanceParam = 0.2)
```

#### **Arguments**

file The input data. More specifically, a data frame with two dimensions, each di-

mension representing the intensity for one color channel.

sensitivity A number between 0.1 and 2 determining sensitivity of the initial clustering,

e.g. the number of clusters. A higher value means more clusters are being

found. Standard is 1.

numOfMarkers The number of primary clusters that are expected according the experiment set

up.

missingClusters

A vector containing the number of primary clusters, which are missing in this

dataset according to the template.

similarityParam

If the distance of a droplet between two or more clusters is very similar, it will not be counted for either. The standard it 0.95, i.e. at least 95% similarity. A sensible value lies between 0 and 1, where 0 means none of the 'rain' droplets

will be counted and 1 means all droplets will be counted.

distanceParam When assigning rain between two clusters, typically the bottom 20% are as-

signed to the lower cluster and the remaining 80% to the higher cluster. This parameter changes the ratio, i.e. a value of 0.1 would assign only 10% to the

lower cluster.

runPeaks 11

#### Value

data The original input data minus the removed events (for plotting)

counts The droplet count for each cluster. firstClusters The position of the primary clusters.

partition The cluster numbers as a CLUE partition (see clue package for more informa-

tion).

# **Examples**

runPeaks

Find the clusters using flowPeaks

#### **Description**

Find the rain and assign it based on the distance to vector lines connecting the cluster centres.

#### Usage

```
runPeaks(file, sensitivity = 1, numOfMarkers, missingClusters = NULL,
    similarityParam = 0.95, distanceParam = 0.2)
```

#### **Arguments**

file The input data. More specifically, a data frame with two dimensions, each di-

mension representing the intensity for one color channel.

sensitivity A number between 0.1 and 2 determining sensitivity of the initial clustering,

e.g. the number of clusters. A higher value means more clusters are being

found. Standard is 1.

numOfMarkers The number of primary clusters that are expected according the experiment set

up.

missingClusters

A vector containing the number of primary clusters, which are missing in this dataset according to the template.

similarityParam

If the distance of a droplet between two or more clusters is very similar, it will not be counted for either. The standard it 0.95, i.e. at least 95% similarity. A sensible value lies between 0 and 1, where 0 means none of the 'rain' droplets will be counted and 1 means all droplets will be counted.

12 runSam

distanceParam When assigning rain between two clusters, typically the bottom 20% are as-

signed to the lower cluster and the remaining 80% to the higher cluster. This parameter changes the ratio, i.e. a value of 0.1 would assign only 10% to the

lower cluster.

# Value

data The original input data minus the removed events (for plotting)

counts The droplet count for each cluster. firstClusters The position of the primary clusters.

partition The cluster numbers as a CLUE partition (see clue package for more informa-

tion).

#### **Examples**

runSam

Find the clusters using SamSPECTRAL

# **Description**

Find the rain and assign it based on the distance to vector lines connecting the cluster centres.

# Usage

```
runSam(file, sensitivity = 1, numOfMarkers, missingClusters = NULL,
  similarityParam = 0.95, distanceParam = 0.2)
```

#### **Arguments**

file The input data. More specifically, a data frame with two dimensions, each di-

mension representing the intensity for one color channel.

sensitivity A number between 0.1 and 2 determining sensitivity of the initial clustering,

e.g. the number of clusters. A higher value means more clusters are being

found. Standard is 1.

numOfMarkers The number of primary clusters that are expected according the experiment set

up.

shearCorrection 13

missingClusters

A vector containing the number of primary clusters, which are missing in this dataset according to the template.

similarityParam

If the distance of a droplet between two or more clusters is very similar, it will not be counted for either. The standard it 0.95, i.e. at least 95% similarity. A sensible value lies between 0 and 1, where 0 means none of the 'rain' droplets will be counted and 1 means all droplets will be counted.

distanceParam

When assigning rain between two clusters, typically the bottom 20% are assigned to the lower cluster and the remaining 80% to the higher cluster. This parameter changes the ratio, i.e. a value of 0.1 would assign only 10% to the lower cluster.

#### Value

data The original input data minus the removed events (for plotting)

counts The droplet count for each cluster. firstClusters The position of the primary clusters.

partition The cluster numbers as a CLUE partition (see clue package for more informa-

tion).

# **Examples**

shearCorrection

Correct for DNA shearing

# **Description**

Longer DNA templates produce a lower droplet count due to DNA shearing. This function normalizes the ddPCRclust result based on a stable marker of different lengths to negate the effect of differences in the lengths of the actual markers of interest. (Work in progress)

# Usage

```
shearCorrection(counts, lengthControl, stableControl)
```

shearCorrection

# Arguments

counts The counts per marker as provided by calculateCPDs.

lengthControl The name of the length Control. If the template name is for example CPT2, the

name in the template should be CPT2-125, where 125 represents the number of

basepairs.

stableControl The name of the stable Control used as a reference for this experiment.

# Value

A linear regression model fitting the length vs ln(ratio) (see lm for details on linear regression).

# **Index**

```
calculateCPDs, 2, 5, 14
cl_{medoid}, 3
createEnsemble, 3, 5
ddPCRclust, 4, 5-8
ddPCRclust-package (ddPCRclust), 4
exportPlots, 5, 6
exportToCSV, 5, 7
exportToExcel, 5, 8
ggsave, 6
lm, 14
parallel, 5
readFiles, 4, 9
readTemplate, 4, 9
runDensity, 5, 10
runPeaks, 5, 11
runSam, 5, 12
shearCorrection, 13
```